

Water Resistant Stopwatch/Clock



Digital Stopwatch

Plus Day/Month/Date Calendar and Alarm

Features:

- Sporty colorful yellow stopwatch with large LCD
- Stopwatch counts up to 23 hours 59 min, 59 seconds accurate to ± 3 seconds/day
- Stopwatch/Chronograph mode with 1/100sec resolution
- Accumulated elapsed time and split time measurements
- Accurately times a two person competition
- Calendar displays day, month, and date
- 12 or 24 hour clock operation with settable hour chime
- Wake and snooze alarm
- Dimensions: 2.3 x 2.8 x 0.5" (57 x 70 x 15mm); Weight: 0.11 lbs (50g)
- Water resistant housing is complete with 39" (1m) neckstrap and battery



Ordering Information:

365510Digital Stopwatch/Clock
365510-NISTDigital Stopwatch/Clock with Calibration Traceable to NIST

